


PABPC3 Antibody (F53970)

Catalog No.	Formulation	Size
F53970-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q9H361
Applications	Western Blot : 1:1000-2000
Limitations	This PABPC3 antibody is available for research use only.



Western blot testing of human 1) Daudi, 2) HepG2, 3) Jurkat, 4) K562, 5) MOLT-4 and 5) Raji cell lysate with PABPC3 antibody at 1:2000. Predicted molecular weight: 70 kDa.

Description

Polyadenylate binding protein 3 binds the poly(A) tail of mRNA. May be involved in cytoplasmic regulatory processes of mRNA metabolism. Binds poly(A) with a slightly lower affinity as compared to PABPC1.

Application Notes

Titration of the PABPC3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 84-117 from the human protein was used as the immunogen for the PABPC3 antibody.

Storage

Aliquot the PABPC3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.