


## PACAP Antibody (R34103)

Catalog No.	Formulation	Size
R34103-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Predicted Reactivity</b>	Cow, Dog, Mouse, Pig, Rat
<b>Format</b>	Antigen affinity purified
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	116
<b>Applications</b>	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:2000
<b>Limitations</b>	This PACAP antibody is available for research use only.



Western blot testing of human hippocampus lysate with PACAP antibody at 1ug/ml.  
Predicted molecular weight: ~19kDa, but can be observed at ~25kDa.

### Description

Additional name(s) for this target protein: Adenylate cyclase activating polypeptide 1; ADCYAP1

### Application Notes

Optimal dilution of the PACAP antibody should be determined by the researcher.

### Immunogen

Amino acids GKRYKQRVKNKGRR were used as the immunogen for this PACAP antibody.

## **Storage**

Aliquot and store the PACAP antibody at -20oC.