


PACAP Antibody (R34103)

Catalog No.	Formulation	Size
R34103-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	116
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This PACAP antibody is available for research use only.



Western blot testing of human hippocampus lysate with PACAP antibody at 1ug/ml.
Predicted molecular weight: ~19kDa, but can be observed at ~25kDa.

Description

Additional name(s) for this target protein: Adenylate cyclase activating polypeptide 1; ADCYAP1

Application Notes

Optimal dilution of the PACAP antibody should be determined by the researcher.

Immunogen

Amino acids GKRYKQRVKNKGRR were used as the immunogen for this PACAP antibody.

Storage

Aliquot and store the PACAP antibody at -20oC.