


PACSIN3 Antibody (R33153)

Catalog No.	Formulation	Size
R33153-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	29763
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This PACSIN3 antibody is available for research use only.



Western blot testing of H460 lysate with PACSIN3 antibody at 1ug/ml. Predicted molecular weight: ~48kDa.

Description

Additional name(s) for this target protein: Protein kinase C and casein kinase substrate in neurons 3

Application Notes

Optimal dilution of the PACSIN3 antibody should be determined by the researcher.

Immunogen

Amino acids APEEDAGGEALGGS were used as the immunogen for this PACSIN3 antibody.

Storage

Aliquot and store the PACSIN3 antibody at -20oC.