


## PAF1 Antibody (R34768)

Catalog No.	Formulation	Size
R34768-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5828
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This PAF1 antibody is available for research use only.



Western blot testing of HeLa lysate with PAF1 antibody at 0.5ug/ml. Predicted/observed molecular weight: ~35kDa.

## Description

Additional name(s) for this target protein: Peroxisome assembly factor-1, peroxisomal membrane protein 3, peroxin-2; PXMP3

## Application Notes

Optimal dilution of the PAF1 antibody should be determined by the researcher.

## **Immunogen**

Amino acids NATVGQSVLNKYKN were used as the immunogen for this PAF1 antibody.

## **Storage**

Aliquot and store the PAF1 antibody at -20oC.