


Pan TRK Antibody / TrkA/B/C [clone RM423] (R20437)

Catalog No.	Formulation	Size
R20437-0.1ML	Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM423
Purity	Protein A purified from animal origin-free supernatant
UniProt	P04629, Q16620, Q16288
Applications	Immunohistochemistry (FFPE) : 1:50-1:100 Western Blot : 1:1000-1:2000
Limitations	This recombinant Pan TRK antibody is available for research use only.



IHC staining of FFPE human brain tissue with recombinant Pan TRK antibody at 1:100.



Western blot testing of mouse brain tissue lysate with recombinant Pan TRK antibody at 1:2000 dilution. Expected molecular weight: 85~140 kDa depending on glycosylation level.

Description

This antibody reacts to human TrkA, TrkB and TrkC.

Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant Pan TRK antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A peptide corresponding to the C-terminus of human TRK was used as the immunogen for the Pan TRK antibody.

Storage

Store the recombinant Pan TRK antibody at -20oC.