


Pard3b Antibody (R35698)

Catalog No.	Formulation	Size
R35698-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	72823
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Pard3b antibody is available for research use only.



Western blot testing of mouse fetal kidney lysate with Pard3b antibody at 0.5ug/ml. Predicted molecular weight: ~133kDa. Both observed bands are blocked by addition of the immunizing peptide.

Description

Additional name(s) for this target protein: Par-3 partitioning defective 3 homolog B

Application Notes

Optimal dilution of the Pard3b antibody should be determined by the researcher.

Immunogen

Amino acids RNQDPRQKNP were used as the immunogen for this Pard3b antibody.

Storage

Aliquot and store the Pard3b antibody at -20oC.