


Pard3b Antibody (R35698)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R35698-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|----------------------|--|
| Availability | 1-3 business days |
| Species Reactivity | Mouse |
| Predicted Reactivity | Rat |
| Format | Antigen affinity purified |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 72823 |
| Applications | Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000 |
| Limitations | This Pard3b antibody is available for research use only. |



Western blot testing of mouse fetal kidney lysate with Pard3b antibody at 0.5ug/ml. Predicted molecular weight: ~133kDa. Both observed bands are blocked by addition of the immunizing peptide.

Description

Additional name(s) for this target protein: Par-3 partitioning defective 3 homolog B

Application Notes

Optimal dilution of the Pard3b antibody should be determined by the researcher.

Immunogen

Amino acids RNQDPRQKNP were used as the immunogen for this Pard3b antibody.

Storage

Aliquot and store the Pard3b antibody at -20oC.