


PATJ Antibody (R34996)

Catalog No.	Formulation	Size
R34996-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10207
Applications	Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This PATJ antibody is available for research use only.



IHC testing of FFPE human testis with PATJ antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining. Spermatid staining is seen.

Description

Additional name(s) for this target protein: PALS1-associated tight junction protein; INADL

Application Notes

Optimal dilution of the PATJ antibody should be determined by the researcher.

Immunogen

Amino acids PENPATDKLQVLQ were used as the immunogen for this PATJ antibody.

Storage

Aliquot and store the PATJ antibody at -20oC.