


PAUF Antibody (F40007)

Catalog No.	Formulation	Size
F40007-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40007-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	C3PTT6
Applications	Western Blot : 1:1000
Limitations	This PAUF antibody is available for research use only.



Description

Pancreatic adenocarcinoma upregulated factor

Application Notes

Titration of the PAUF antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 109-137 from the human protein was used as the immunogen for this PAUF antibody.

Storage

Aliquot the PAUF antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.