


PAX3 Antibody (Biotin Conjugate) (R36304BTN)

Catalog No.	Formulation	Size
R36304BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Rat
Format	Biotin Conjugate
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	P23760
Gene ID	5077
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This PAX3 antibody is available for research use only.



Western blot testing of human duodenum lysate with biotinylated PAX3 antibody at 0.5ug/ml. Predicted molecular weight: can detect 7 isoforms ranging from 22~56kDa.

Description

Additional name(s) for this target protein: Paired box 3

Application Notes

Optimal dilution of the PAX3 antibody should be determined by the researcher.

1. This Ab can detect isoforms 1-7, ranging from 206-505 amino acids in length.

Immunogen

Amino acids TTLAGAVPRMM were used as the immunogen for this PAX3 antibody. (1)

Storage

Aliquot and store the PAX3 antibody at -20oC.