


PAX3 Antibody (R36091)

Catalog No.	Formulation	Size
R36091-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5077
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This PAX3 antibody is available for research use only.



Western blot testing of mouse brain lysate with PAX3 antibody at 2ug/ml. Predicted molecular weight: can detect 6 isoforms ranging from 44~56kDa.

Description

Additional name(s) for this target protein: Paired box 3

Application Notes

Optimal dilution of the PAX3 antibody should be determined by the researcher.

1. This Ab can detect isoforms ranging from 403-505 amino acids in length.

Immunogen

Amino acids STVHQSTIPSNPD were used as the immunogen for this PAX3 antibody.

Storage

Aliquot and store the PAX3 antibody at -20oC.