


PAX5 Antibody (R33523)


Catalog No.	Formulation	Size
R33523-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5079
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This PAX5 antibody is available for research use only.



Western blot testing of human lymph lysate with PAX5 antibody at 0.3ug/ml. Predicted molecular weight ~42 kDa.



Western blot testing of 1) mouse spleen lysate at 1ug/ml and 2) rat spleen lysate at 0.3ug/ml with PAX5 antibody. Predicted molecular weight ~42 kDa.

Description

Additional name(s) for this target protein: Paired box gene 5

Application Notes

Optimal dilution of the PAX5 antibody should be determined by the researcher.

Immunogen

Amino acids DLEKNYPTPRTSR were used as the immunogen for this PAX5 antibody.

Storage

Aliquot and store the PAX5 antibody at -20oC.