

PAX6 Antibody / Paired box protein Pax-6 [clone MSVA-706M] (V6102)

Catalog No.	Formulation	Size
V6102-100UG	Antibody in 1X PBS with 0.05% BSA, 0.05% sodium azide	100 ug
V6102-20UG	Antibody in 1X PBS with 0.05% BSA, 0.05% sodium azide	20 ug

[Bulk quote request](#)

Species Reactivity	Human
Format	Purified
Host	Mouse
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG1, kappa
Clone Name	MSVA-706M
UniProt	P26367
Localization	Nucleus
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml
Limitations	This PAX6/Paired box protein Pax-6 antibody is available for research use only.

Application Notes

1. Optimal dilution of the PAX6/Paired box protein Pax-6 antibody should be determined by the researcher.
2. This PAX6/Paired box protein Pax-6 antibody is recombinantly produced by expression in CHO cells.

Immunogen

A recombinant fragment (N-terminus; aa 1-300) of human PAX6 protein was used as the immunogen for the PAX6/Paired box protein Pax-6 antibody.

Storage

PAX6/Paired box protein Pax-6 antibody with sodium azide - store at 2 to 8°C; antibody without sodium azide - store at -20 to -80°C.

