

PAX8 Antibody / Paired box protein Pax-8 [clone PAX8/12570] (V6019)

Catalog No.	Formulation	Size
V6019-100UG	0.2 mg/ml in 1X PBS with 0.05% BSA, 0.05% sodium azide	100 ug
V6019-20UG	0.2 mg/ml in 1X PBS with 0.05% BSA, 0.05% sodium azide	20 ug
V6019SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

[Bulk quote request](#)

Species Reactivity	Human
Format	Purified
Host	Mouse
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG1, kappa
Clone Name	PAX8/12570
UniProt	Q06710
Localization	Nucleus
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml
Limitations	This PAX8/Paired box protein Pax-8 antibody is available for research use only.

Application Notes

Optimal dilution of the PAX8/Paired box protein Pax-8 antibody should be determined by the researcher.

Immunogen

A recombinant fragment (around amino acids 50-300) of human PAX8 protein (exact sequence is proprietary) was used as the immunogen for the PAX8/Paired box protein Pax-8 antibody.

Storage

PAX8/Paired box protein Pax-8 antibody with sodium azide - store at 2 to 8°C; antibody without sodium azide - store at -20 to -80°C.

