


## Pbr Antibody (R35566)

Catalog No.	Formulation	Size
R35566-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	12257
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This Pbr antibody is available for research use only.



Western blot testing of mouse kidney lysate [A] and of mouse spleen lysate [B] with Pbr antibody at 1ug/ml. Predicted molecular weight: ~19 kDa.

## Description

Additional name(s) for this target protein: Peripheral benzodiazepine receptor, translocator protein (18kDa), Tspo

## Application Notes

Optimal dilution of the Pbr antibody should be determined by the researcher.

## Immunogen

Amino acids RDNSGRRGGSRLPE were used as the immunogen for this Pbr antibody.

## **Storage**

Aliquot and store the Pbr antibody at -20oC.