


Pbr Antibody (R35566)

Catalog No.	Formulation	Size
R35566-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	12257
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This Pbr antibody is available for research use only.



Western blot testing of mouse kidney lysate [A] and of mouse spleen lysate [B] with Pbr antibody at 1ug/ml. Predicted molecular weight: ~19 kDa.

Description

Additional name(s) for this target protein: Peripheral benzodiazepine receptor, translocator protein (18kDa), Tspo

Application Notes

Optimal dilution of the Pbr antibody should be determined by the researcher.

Immunogen

Amino acids RDNSGRRGGSRLPE were used as the immunogen for this Pbr antibody.

Storage

Aliquot and store the Pbr antibody at -20oC.