


PCBP1 Antibody (R33113)

Catalog No.	Formulation	Size
R33113-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5093
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This PCBP1 antibody is available for research use only.



Western blot testing of HeLa lysate with PCBP1 antibody at 0.01ug/ml. Predicted molecular weight: ~38kDa.

Description

Additional name(s) for this target protein: Poly(rC) binding protein 1

Application Notes

Optimal dilution of the PCBP1 antibody should be determined by the researcher.

Immunogen

Amino acids ADLAKLNQVARQQSH were used as the immunogen for this PCBP1 antibody.

Storage

Aliquot and store the PCBP1 antibody at -20oC.