


## PCBP4 Antibody (R34810)


Catalog No.	Formulation	Size
R34810-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	57060
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This PCBP4 antibody is available for research use only.



Western blot testing of HeLa cell lysate with PCBP4 antibody at 0.03ug/ml. Predicted molecular weight: ~41kDa.



Western blot of HEK293 lysate overexpressing PCBP4 probed with PCBP4 antibody (mock transfection in lane 1). Predicted molecular weight: ~41kDa.

## Description

Additional name(s) for this target protein: Poly(rC) binding protein 4

## Application Notes

Optimal dilution of the PCBP4 antibody should be determined by the researcher.

## Immunogen

Amino acids NGSKKAERQKFSPY were used as the immunogen for this PCBP4 antibody.

## Storage

Aliquot and store the PCBP4 antibody at -20oC.