


PCF11 Antibody (RQ8478)

Catalog No.	Formulation	Size
RQ8478	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

[Bulk quote request](#)

Availability	1-3 days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	O94913
Applications	Western Blot : 0.5-1ug/ml ELISA : 0.1-0.5ug/ml
Limitations	This PCF11 antibody is available for research use only.



Western blot testing of 1) human K562, 2) human MCF7, 3) rat testis, 4) rat C6 and 5) mouse testis tissue lysate with PCF11 antibody. Predicted molecular weight ~173 kDa.

Description

The protein encoded by this gene binds to CLP1 to form pre-mRNA cleavage factor II_m. The encoded protein is necessary for efficient Pol II transcription termination and may be involved in degradation of the 3' product of polyA site cleavage.

Application Notes

Optimal dilution of the PCF11 antibody should be determined by the researcher.

Immunogen

An E.coli-derived human recombinant protein (amino acids N1214-V1555) was used as the immunogen for the PCF11 antibody.

Storage

After reconstitution, the PCF11 antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.