


PCGF2 Antibody (R34534)

Catalog No.	Formulation	Size
R34534-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7703
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This PCGF2 antibody is available for research use only.



Western blot testing of human lung lysate with PCGF2 antibody at 1ug/ml. Predicted molecular weight: ~38kDa.

Description

Additional name(s) for this target protein: Polycomb group ring finger 2, MEL18

Application Notes

Optimal dilution of the PCGF2 antibody should be determined by the researcher.

Immunogen

Amino acids KMTVNGAPVPPLT were used as the immunogen for this PCGF2 antibody.

Storage

Aliquot and store the PCGF2 antibody at -20oC.