


PCGF3 Antibody (R35058)

Catalog No.	Formulation	Size
R35058-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10336
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This PCGF3 antibody is available for research use only.



Western blot testing of human heart lysate and PCGF3 antibody at 1ug/ml. Predicted molecular weight: ~27kDa.

Description

Additional name(s) for this target protein: Ring finger protein 3, RNF3, polycomb group ring finger 3

Application Notes

Optimal dilution of the PCGF3 antibody should be determined by the researcher.

1. This PCGF3 antibody will detect the 247 AA isoform. It will not detect the 132 AA isoform.

Immunogen

Amino acids PLLLHYRPMKMDLL were used as the immunogen for this PCGF3 antibody.

Storage

Aliquot and store the PCGF3 antibody at -20°C.