


## PCNA Antibody Biotin Conjugate (R34806BTN)

Catalog No.	Formulation	Size
R34806BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human, Mouse, Rat, Pig
<b>Predicted Reactivity</b>	Dog, Cow
<b>Format</b>	Biotin Conjugate
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	5111
<b>Applications</b>	Western Blot : 1-3ug/ml Enzyme ImmunoAssay : binds recombinant protein ELISA (peptide) LOD : 1:128000
<b>Limitations</b>	This PCNA antibody is available for research use only.



Western blot testing of HeLa lysate with biotinylated PCNA antibody at 1ug/ml. Expected molecular weight: 29-36 kDa.

### Description

Additional name(s) for this target protein: Proliferating cell nuclear antigen

### Application Notes

Optimal dilution of the PCNA antibody should be determined by the researcher.

## **Immunogen**

Amino acids NGNIKLSQTSNVD were used as the immunogen for this PCNA antibody.

## **Storage**

Aliquot and store the PCNA antibody at -20oC.