

## PCNA Antibody Biotin Conjugate (R34806BTN)

Catalog No.	Formulation	Size
R34806BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat, Pig
Predicted Reactivity	Dog, Cow
Format	Biotin Conjugate
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5111
Applications	Western Blot : 1-3ug/ml Enzyme ImmunoAssay : binds recombinant protein ELISA (peptide) LOD : 1:128000
Limitations	This PCNA antibody is available for research use only.



Western blot testing of HeLa lysate with biotinylated PCNA antibody at 1ug/ml. Expected molecular weight: 29-36 kDa.

### Description

Additional name(s) for this target protein: Proliferating cell nuclear antigen

### Application Notes

Optimal dilution of the PCNA antibody should be determined by the researcher.

## **Immunogen**

Amino acids NGNIKLSQTSNVD were used as the immunogen for this PCNA antibody.

## **Storage**

Aliquot and store the PCNA antibody at -20oC.