


PCNXL2 Antibody (F53858)

Catalog No.	Formulation	Size
F53858-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F53858-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	A6NKB5
Applications	Western Blot : 1:1000-1:2000
Limitations	This PCNXL2 antibody is available for research use only.



Western blot testing of human 1) HeLa, 2) HT-29 and 3) WiDr cell lysate with PCNXL2 antibody at 1:2000. Predicted molecular weight: 237 kDa.

Description

Pecanex-like protein 2 may play a role in tumorigenesis of colorectal carcinomas with high microsatellite instability (MSI-H).

Application Notes

Titration of the PCNXL2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 467-499 from the human protein was used as the immunogen for the PCNXL2 antibody.

Storage

Aliquot the PCNXL2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.