


PCSK6 Antibody / PACE4 (R36373)

Catalog No.	Formulation	Size
R36373-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5046
Applications	Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This PCSK6 antibody is available for research use only.



IHC testing of FFPE human skeletal muscle and PCSK6 antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Proprotein convertase subtilisin/kexin type 6

Application Notes

Optimal dilution of the PCSK6 antibody should be determined by the researcher.

1. This PCSK6 antibody is specific for isoforms a/b/g/h.

Immunogen

Amino acids YHTFSAHQSRSR were used as the immunogen for this PCSK6 antibody.

Storage

Aliquot and store the PCSK6 antibody at -20oC.