


## PCSK9 Antibody (F52952)

Catalog No.	Formulation	Size
F52952-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F52952-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

### Bulk quote request

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Antigen affinity
<b>UniProt</b>	Q8NBP7
<b>Applications</b>	Western Blot : 1:1000
<b>Limitations</b>	This PCSK9 antibody is available for research use only.



Western blot analysis of lysate from SK-N-MC cell line using PCSK9 antibody at 1:1000.  
Predicted size: Pro/mature ~74/64 kDa

## Description

Crucial player in the regulation of plasma cholesterol homeostasis. Binds to low-density lipoprotein receptor family members: low density lipoprotein receptor (LDLR), very low density lipoprotein receptor (VLDLR), apolipoprotein E receptor (LRP1/APOER) and apolipoprotein receptor 2 (LRP8/APOER2), and promotes their degradation in intracellular acidic compartments. Acts via a non-proteolytic mechanism to enhance the degradation of the hepatic LDLR through a clathrin LDLRAP1/ARH-mediated pathway. May prevent the recycling of LDLR from endosomes to the cell surface or direct it to lysosomes for degradation. Can induce ubiquitination of LDLR leading to its subsequent degradation. Inhibits intracellular degradation of APOB via the autophagosome/lysosome pathway in a LDLR-independent manner. Involved in the disposal

of non-acetylated intermediates of BACE1 in the early secretory pathway. Inhibits epithelial Na(+) channel (ENaC)-mediated Na(+) absorption by reducing ENaC surface expression primarily by increasing its proteasomal degradation. Regulates neuronal apoptosis via modulation of LRP8/APOER2 levels and related anti-apoptotic signaling pathways.

## Application Notes

Titration of the PCSK9 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

This PCSK9 antibody was produced from a rabbit immunized with a KLH conjugated synthetic peptide between 654-688 amino acids from the C-terminal region of human PCSK9.

## Storage

Aliquot the PCSK9 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.