


PDCD4 Antibody (R33777)


Catalog No.	Formulation	Size
R33777-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Pig, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	27250
Localization	Nuclear, cytoplasmic
Applications	Western Blot : 0.1-0.3ug/ml IHC (FFPE) : 5ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This PDCD4 antibody is available for research use only.



Western blot testing of HeLa lysate with PDCD4 antibody at 0.1ug/ml. Expected molecular weight: 50-60 kDa.



IHC testing of FFPE human prostate tissue with PDCD4 antibody at 5ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

Description

Additional name(s) for this target protein: Programmed cell death 4

Application Notes

Optimal dilution of the PDCD4 antibody should be determined by the researcher.

Immunogen

Amino acids EGKASHREMTSK were used as the immunogen for this PDCD4 antibody.

Storage

Aliquot and store the PDCD4 antibody at -20oC.