


## Pdgfra Antibody (F43668)

Catalog No.	Formulation	Size
F43668-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F43668-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Mouse
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Antigen affinity
<b>UniProt</b>	P26618
<b>Applications</b>	Western Blot : 1:1000
<b>Limitations</b>	This Pdgfra antibody is available for research use only.



### Description

Receptor that binds both PDGFA and PDGFB and has a tyrosine-protein kinase activity.

### Application Notes

Titration of the Pdgfra antibody may be required due to differences in protocols and secondary/substrate sensitivity.

### Immunogen

A portion of amino acids 678-706 from the mouse protein was used as the immunogen for this Pdgfra antibody.

## **Storage**

Aliquot the Pdgfra antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.