


Pdgfra Antibody (F43668)

Catalog No.	Formulation	Size
F43668-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F43668-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P26618
Applications	Western Blot : 1:1000
Limitations	This Pdgfra antibody is available for research use only.



Description

Receptor that binds both PDGFA and PDGFB and has a tyrosine-protein kinase activity.

Application Notes

Titration of the Pdgfra antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 678-706 from the mouse protein was used as the immunogen for this Pdgfra antibody.

Storage

Aliquot the Pdgfra antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.