


PDIA2 Antibody Biotin Conjugate (R34048BTN)

Catalog No.	Formulation	Size
R34048BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Biotin Conjugate
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	64714
Applications	Western Blot : 0.03-0.1ug/ml Immunohistochemistry (FFPE) : suitable ELISA (peptide) LOD : 1:64000
Limitations	This PDIA2 antibody is available for research use only.



Western blot testing of rat pancreas lysate with biotinylated PDIA2 antibody at 0.03ug/ml. Predicted molecular weight ~58/70kDa (unmodified/glycosylated).

Description

Additional name(s) for this target protein: Protein disulfide isomerase family A, member 2

Application Notes

Optimal dilution of the PDIA2 antibody should be determined by the researcher.

Immunogen

Amino acids EYKSTRDLETFSK were used as the immunogen for this PDIA2 antibody.

Storage

Aliquot and store the PDIA2 antibody at -20oC.