


PDIA2 Antibody (R34048)


Catalog No.	Formulation	Size
R34048-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	64714
Applications	Western Blot : 0.01-0.03ug/ml Immunohistochemistry (FFPE) : 3-5ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This PDIA2 antibody is available for research use only.



Western blot testing of human [A], mouse [B] and rat [C] pancreas lysates with PDIA2 antibody at 0.1ug/ml. Predicted molecular weight ~58/70kDa (unmodified/glycosylated).



IHC staining of FFPE human pancreas with PDIA2 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Protein disulfide isomerase family A, member 2

Application Notes

Optimal dilution of the PDIA2 antibody should be determined by the researcher.

Immunogen

Amino acids EYKSTRDLETFSK were used as the immunogen for this PDIA2 antibody.

Storage

Aliquot and store the PDIA2 antibody at -20oC.