


## PDLIM2 Antibody (R34886)

Catalog No.	Formulation	Size
R34886-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	64236
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This PDLIM2 antibody is available for research use only.



Western blot testing of mouse lung lysate with PDLIM2 antibody at 0.3ug/ml.  
Predicted/observed molecular weight: ~38kDa.

## Description

Additional name(s) for this target protein: PDZ and LIM domain 2, MYSTIQUE, SLIM

## Application Notes

Optimal dilution of the PDLIM2 antibody should be determined by the researcher.

1. This PDLIM2 antibody will detect isoforms 1 and 2.

## **Immunogen**

Amino acids NQAVRIQEGRYRHP were used as the immunogen for this PDLIM2 antibody.

## **Storage**

Aliquot and store the PDLIM2 antibody at -20oC.