


PDLIM4 Antibody (R33674)


Catalog No.	Formulation	Size
R33674-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8572
Localization	Cytoplasmic
Applications	Western Blot : 0.1-0.5ug/ml IHC (FFPE) : 3-5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This PDLIM4 antibody is available for research use only.



Western blot testing of A431 lysate with PDLIM4 antibody at 0.2ug/ml. Predicted molecular weight: ~35 kDa.



IHC testing of FFPE human breast tissue with PDLIM4 antibody. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

Description

Additional name(s) for this target protein: PDZ and LIM domain 4

Application Notes

Optimal dilution of the PDLIM4 antibody should be determined by the researcher.

Immunogen

Amino acids DVVAVYPNAKVELV were used as the immunogen for this PDLIM4 antibody.

Storage

Aliquot and store the PDLIM4 antibody at -20oC.