


PDXP Antibody (R35288)


Catalog No.	Formulation	Size
R35288-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	57026
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This PDXP antibody is available for research use only.



Western blot testing of human heart lysate with PDXP antibody at 0.1ug/ml. Predicted molecular weight: ~32kDa.



Western blot of mouse (A) and rat (B) brain lysates with PDXP antibody at 0.1ug/ml. Predicted molecular weight: ~32kDa.

Description

Additional name(s) for this target protein: Pyridoxal phosphatase; Chronophin

Application Notes

Optimal dilution of the PDXP antibody should be determined by the researcher.

Immunogen

Amino acids QHDLVPHYVES were used as the immunogen for this PDXP antibody.

Storage

Aliquot and store the PDXP antibody at -20oC.