

Periostin Antibody (R36408)

Catalog No.	Formulation	Size
R36408-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10631
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Periostin antibody is available for research use only.



Western blot testing of human breast cancer lysate with Periostin antibody at 0.1ug/ml.
Predicted molecular weight: ~93kDa, observed here at ~85kDa.

Description

Additional name(s) for this target protein: POSTN; Osteoblast-specific factor 2 (OSF2)

Application Notes

Optimal dilution of the Periostin antibody should be determined by the researcher.

Immunogen

Amino acids TTQRYSDASKLREE were used as the immunogen for this Periostin antibody.

Storage

Aliquot and store the Periostin antibody at -20oC.