

Peroxiredoxin 2 Antibody (R34939)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R34939-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|----------------------|---|
| Availability | 1-3 business days |
| Species Reactivity | Human, Mouse |
| Predicted Reactivity | Rat |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 7001 |
| Applications | Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:16000 |
| Limitations | This Peroxiredoxin 2 antibody is available for research use only. |



Western blot testing of human brain lysate with Peroxiredoxin 2 antibody at 1ug/ml.
Predicted molecular weight: ~22kDa.

Description

Additional name(s) for this target protein: PRDX2

Application Notes

Optimal dilution of the Peroxiredoxin 2 antibody should be determined by the researcher.

1. This Peroxiredoxin 2 antibody will detect isoforms a and b.

Immunogen

Amino acids NVDDSKEYFSKHN were used as the immunogen for this Peroxiredoxin 2 antibody.

Storage

Aliquot and store the Peroxiredoxin 2 antibody at -20oC.