

## Peroxiredoxin 2 Antibody (R34939)

| Catalog No.  | Formulation   | Size   |
|--------------|---|--------|
| R34939-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                      |   |
|----------------------|---|
| Availability         | 1-3 business days   |
| Species Reactivity   | Human, Mouse  |
| Predicted Reactivity | Rat   |
| Format               | Antigen affinity purified   |
| Clonality            | Polyclonal (goat origin)  |
| Isotype              | Goat Ig   |
| Purity               | Antigen affinity  |
| Gene ID              | 7001  |
| Applications         | Western Blot : 0.3-1ug/ml<br>ELISA (peptide) LOD : 1:16000        |
| Limitations          | This Peroxiredoxin 2 antibody is available for research use only. |



Western blot testing of human brain lysate with Peroxiredoxin 2 antibody at 1ug/ml.  
Predicted molecular weight: ~22kDa.

## Description

Additional name(s) for this target protein: PRDX2

## Application Notes

## Immunogen

Amino acids NVDDSKEYFSKHN were used as the immunogen for this Peroxiredoxin 2 antibody.

## Storage

Aliquot and store the Peroxiredoxin 2 antibody at -20oC.