


Peroxiredoxin 6 Antibody (R36124)

Catalog No.	Formulation	Size
R36124-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Pig
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9588
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Peroxiredoxin 6 antibody is available for research use only.



Western blot testing of human heart [A] and liver [B] lysates with Peroxiredoxin 6 antibody at 0.3ug/ml. Expected molecular weight: ~25kDa.



Western blot of pig liver lysate with Peroxiredoxin 6 antibody at 0.1ug/ml. Expected molecular weight: ~25kDa.

Description

Additional name(s) for this target protein: PRDX6

Application Notes

Optimal dilution of the Peroxiredoxin 6 antibody should be determined by the researcher.

Immunogen

Amino acids TAEKRVATPVD were used as the immunogen for this Peroxiredoxin 6 antibody.

Storage

Aliquot and store the Peroxiredoxin 6 antibody at -20oC.