


## PERP Antibody (R33273)


Catalog No.	Formulation	Size
R33273-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	64065
Applications	Western Blot : 1-3ug/ml IHC (FFPE) : 5ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This PERP antibody is available for research use only.



Western blot testing of human liver lysate with PERP antibody at 1ug/ml. Predicted molecular weight: ~21 kDa.



IHC testing of FFPE human skin tissue with PERP antibody at 5ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

## Description

Additional name(s) for this target protein: P53 apoptosis effector related to PMP-22

## Application Notes

Optimal dilution of the PERP antibody should be determined by the researcher.

## Immunogen

Amino acids CLPNYEDDLLGNAK were used as the immunogen for this PERP antibody.

## Storage

Aliquot and store the PERP antibody at -20oC.