


PERP Antibody (R33273)


Catalog No.	Formulation	Size
R33273-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	64065
Applications	Western Blot : 1-3ug/ml IHC (FFPE) : 5ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This PERP antibody is available for research use only.



Western blot testing of human liver lysate with PERP antibody at 1ug/ml. Predicted molecular weight: ~21 kDa.



IHC testing of FFPE human skin tissue with PERP antibody at 5ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

Description

Additional name(s) for this target protein: P53 apoptosis effector related to PMP-22

Application Notes

Optimal dilution of the PERP antibody should be determined by the researcher.

Immunogen

Amino acids CLPNYEDDLLGNAK were used as the immunogen for this PERP antibody.

Storage

Aliquot and store the PERP antibody at -20oC.