


## Pescadillo Antibody (R34090)


Catalog No.	Formulation	Size
R34090-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23481
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : 2.5-5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Pescadillo antibody is available for research use only.



Western blot testing of HeLa lysate with Pescadillo antibody at 1ug/ml.  
Predicted/observed molecular weight: ~68kDa.



IHC staining of FFPE human testis with Pescadillo antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

## Description

Additional name(s) for this target protein: PES1; pescadillo 1; PES1

## Application Notes

Optimal dilution of the Pescadillo antibody should be determined by the researcher.

## Immunogen

Amino acids GGLEKKKYERGSA were used as the immunogen for this Pescadillo antibody.

## Storage

Aliquot and store the Pescadillo antibody at -20oC.