


Pex1 Antibody (R35193)

Catalog No.	Formulation	Size
R35193-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	71382
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Pex1 antibody is available for research use only.



Western blot testing of mouse liver lysate with Pex1 antibody at 0.3ug/ml. Predicted molecular weight: ~137kDa, observed here at ~150kDa.

Description

Additional name(s) for this target protein: Peroxin 1, Peroxisomal biogenesis factor 1

Application Notes

Optimal dilution of the Pex1 antibody should be determined by the researcher.

Immunogen

Amino acids QAKESTFPKEGDAH were used as the immunogen for this Pex1 antibody.

Storage

Aliquot and store the Pex1 antibody at -20oC.