


PEX13 Antibody (R33634)

Catalog No.	Formulation	Size
R33634-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5194
Applications	Immunohistochemistry (FFPE) : 2-4ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This PEX13 antibody is available for research use only.



IHC testing of FFPE human liver with PEX13 antibody at 3ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: Peroxisome biogenesis factor 13

Application Notes

Optimal dilution of the PEX13 antibody should be determined by the researcher.

Immunogen

Amino acids DSIGKDGEKQDL were used as the immunogen for this PEX13 antibody.

Storage

Aliquot and store the PEX13 antibody at -20oC.