


PEX26 Antibody (R35334)

Catalog No.	Formulation	Size
R35334-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	55670
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This PEX26 antibody is available for research use only.



Western blot testing of human kidney lysate with PEX26 antibody at 0.5ug/ml. Predicted molecular weight: ~34kDa.

Description

Additional name(s) for this target protein: Peroxisomal biogenesis factor 26

Application Notes

Optimal dilution of the PEX26 antibody should be determined by the researcher.

Immunogen

Amino acids QKPNLEGSVSHK were used as the immunogen for this PEX26 antibody.

Storage

Aliquot and store the PEX26 antibody at -20oC.