


PEX6 Antibody (R33930)

Catalog No.	Formulation	Size
R33930-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5190
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This PEX6 antibody is available for research use only.



Western blot testing of rat liver lysate with PEX6 antibody at 0.3ug/ml. Predicted molecular weight: ~104kDa.

Description

Additional name(s) for this target protein: Peroxisomal biogenesis factor 6

Application Notes

Optimal dilution of the PEX6 antibody should be determined by the researcher.

Immunogen

Amino acids EQELLYKRIQRK were used as the immunogen for this PEX6 antibody.

Storage

Aliquot and store the PEX6 antibody at -20oC.