


PGC1A Antibody (R33457)

Catalog No.	Formulation	Size
R33457-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10891
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This PGC1A antibody is available for research use only.



Western blot testing of human skeletal muscle lysate with PGC1A antibody at 0.1ug/ml.
Predicted molecular weight ~91 kDa.

Description

Additional name(s) for this target protein: Peroxisome proliferator-activated receptor gamma, coactivator 1 alpha; PPARGC1A

Application Notes

Optimal dilution of the PGC1A antibody should be determined by the researcher.

Immunogen

Amino acids DGLFDDSEDESDK were used as the immunogen for this PGC1A antibody.

Storage

Aliquot and store the PGC1A antibody at -20oC.