


PGK1 Antibody / Phosphoglycerate kinase 1 (R35940)

Catalog No.	Formulation	Size
R35940-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5230
Localization	Nuclear, cytoplasmic
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This PGK1 antibody is available for research use only.



Western blot testing of HeLa lysate with PGK1 antibody at 0.1ug/ml. Expected/observed molecular weight ~44kDa.

Description

Additional name(s) for this target protein: Phosphoglycerate kinase 1

Application Notes

Optimal dilution of the PGK1 antibody should be determined by the researcher.

Immunogen

Amino acids SKKYAEAVTRAKQ were used as the immunogen for this PGK1 antibody.

Storage

Aliquot and store the PGK1 antibody at -20oC.