


PGRP-S Antibody (R35013)


Catalog No.	Formulation	Size
R35013-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8993
Applications	Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This PGRP-S antibody is available for research use only.



Western blot testing of Jurkat lysate with PGRP-S antibody at 0.3ug/ml. Predicted molecular weight: ~22/28kDa (unmodified/glycosylated).



IHC staining of FFPE human spleen with PGRP-S antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Peptidoglycan recognition protein 1; PGLYRP1

Application Notes

Optimal dilution of the PGRP-S antibody should be determined by the researcher.

Immunogen

Amino acids PGNQLYHLIQN were used as the immunogen for this PGRP-S antibody.

Storage

Aliquot and store the PGRP-S antibody at -20oC.