


PHF11 Antibody (R34472)

Catalog No.	Formulation	Size
R34472-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	51131
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This PHF11 antibody is available for research use only.



Western blot testing of A) Daudi, B) Jurkat whole and C) Jurkat nuclear lysates with PHF11 antibody at 1ug/ml. Predicted/observed molecular weight: ~37kDa.

Description

Additional name(s) for this target protein: PHD finger protein 11

Application Notes

Optimal dilution of the PHF11 antibody should be determined by the researcher.

Immunogen

Amino acids KIHASQQRWQQLKE were used as the immunogen for this PHF11 antibody.

Storage

Aliquot and store the PHF11 antibody at -20oC.